

Analysing Stable Isotope Ratios in HAAs by GC-c-IRMS.

Abstract

Determination of stable isotope ratio ($^{13}\text{C}/^{12}\text{C}$ and $^{15}\text{N}/^{14}\text{N}$) in HAAs by GC-c-IRMS. This analysis can be performed on faunal material as well as on sediment samples. The preparation is given only for biological samples

This is the version of March 2007 by Bart Veuger.

Keywords: fatty acids, chromatogram, NIOO-CEME .

1. Preparation of biological samples

- Mark 12x32 screw cap vials with sample name/code using glass engraving tool (RE lab).
- Pre-fit caps: Put closed (i.e. not open-top), teflon-lined caps on vials and warm in oven or blockheater. When warm: screw caps tight and directly unscrew (for optimal fit)
- 1 - 5 mg fresh (wet) fauna or freeze dried equivalent (0.1 - 0.5 mg). [or hydrolyse more material and later derivatize only part of hydrolysate (in 2nd 12x32 vial)]
- Preferably destruct/homogenize using micro tissue grinder. Put material directly in 6M HCl (or first in Milli-Q and then add same volume 12M HCl). As much liquid as needed but as little as possible (max 1.5 ml in 12x32 vials).
- Create N₂ headspace by flushing with N₂ (5.0) for few seconds
- Vials in blockheater or oven at 110 °C for 20h (overnight)
- Take tubes from blockheater/oven (after 20 h) and allow to cool down.
- Add internal standard (L-Norleucine (Nle) stock solution, 50 mg in 20 ml MQ).
- Indication: 10 ul Nle per mg dry material hydrolyzed (or first dilute Nle stock 10x and then add 100ul/gdw). Be precise and note added volume.
- Evaporate under N₂ (5.0) flow in block heater at 60°C
- Also evaporate some aliquots of standard solution

2. Derivatation

This is the version of September 2009 by Bart Veuger.

- Prepare P-buffer: For 100 ml (100 samples): 0.36 g KH₂PO₄ [A.0.56 VA] and 0.71 g Na₂HPO₄ (*2 H₂O) [A.0.25-S5] in 100 ml Milli-Q
- Switch on block heater (110°C)
- Take samples from fridge (1 series = 24 samples, including standard mix)
- Acidify isopropanol (IP, also called 2 -propanol) by adding acetyl chloride (AC) [A.0.25-K] 5:1 in test tube on ice . For 24 samples: 10 ml IP + 2 ml AC
- Add 0.5 ml acidified IP per sample (glass Pasteur pipette as tip on precision pipette) and close tight
- Samples in block heater at 110°C for 90 min. After few minutes, check in caps are tight (if not, liquid in vial will bubble and volume will decrease)
- Meanwhile, - Label GC vials , Start heating plate with 2 aluminum blocks (50°C = “3”), Prepare evaporation setup (“Octopus” with new Pasteur pipettes) and take P buffer from fridge.
- After 90 min., take vials from heater and let cool down.
- Attach vials to “Octopus” and evaporate under gentle N₂ (5.0) flow at 50°C (= “3”)
- When dry, add 100-200 µl dichloromethane (DCM) [A.0.42-under fume hood] and evaporate (with glass 1 ml syringe, no need to take vials from evaporation setup)
- When dry, add 150 l DCM (with glass 1 ml syringe) and 50 l pentafluoropropionic anhydride (PFA) [A.0.42-under fume hood] (with glass 250 l syringe) and close tight
- Heat vials in block heater for 10 min at 110°C, check if caps are tight.
- After 10 min: cool down to room temperature

3. Solvent extraction

- Add 0.5 ml chloroform and 1 ml P-buffer (with Pasteur pipette as tip on precision pipette)
- Shake/vortex with “Vortex Genie 2” in Molecular lab (max speed for couple of minutes) or shake manually until lower fraction (chloroform) is clear
- Centrifuge (2000 rpm for 10 min)

- Remove lower fraction (chloroform) using 1 ml glass syringe with stainless steel needle (with blunt tip). Hold vial at 45° angle > insert needle along the vial wall > place tip of needle at lowest point > take up chloroform (leave last bit to avoid picking up contamination) > retract needle, > wipe needle with tissue > transfer chloroform to GC vial. Rinse syringe with some chloroform between samples.
- Evaporate chloroform in GC vials. Preferably leave samples overnight in fume hood. If not completely evaporated the next morning: place vials on the side and move liquid from tip to the wall of the vials (increase evaporation surface). Speeding up evaporation by heating and/or N₂ flow is not recommended since derivatized amino acids are very volatile and therefore easily lost by evaporation (especially Ala !)

4. Analysis

- When samples are completely dry: dissolve in ethyl acetate (EA). Minimum volume = 20-50 l. Preferably, start with one sample in 50 l > analyze to determine proper EA volume > dissolve other samples in proper volume. Ideal samples yield highest peaks (usually Asp and Glu) of 5 to 7 V (on GC-c-IRMS in 13C mode), this generally gives a high enough D-Ala peak without other peaks becoming too high (avoiding column overloading effects).
- For GC-c-IRMS: Cap vials with teflon/butyl lined caps (for GC-FID/MS: open caps with Teflon inlays)
- Store samples in freezer until analysis
- Analyze samples as soon after derivatization as possible
- After analysis, recap with new Teflon/butyl lined caps and store in freezer.

Affiliation:

NIOO-CEME